



# Innovation in Injection Techniques

## Cannulas

**01 Complete elimination of a major risk**  
Complete elimination of the risk of intravascular injection during the injection procedure, the consequence of which is the elimination of the risk of necrosis.

**02 Reduction of the invasive aspect of the injections**


- Reduction of the risk of ecchymosis
- The injection becomes less painful thanks to the reduction of the number of pricks on the face. Because the needle doesn't prick and doesn't cut, it slides under the skin without damaging its structure. Before: one prick per injection. Today: through a single point numerous injections can be done following a spoke wheel pattern.

**03 Better application of the product**  
The injections can be done on all areas of the face and body, even in difficult areas. The needle allows us to always control the level of injection. The application of different filler products is done with greater safety, faster, and is of better quality: harmonization of injections and reduction of risks of indentation.



**04 No social exclusion following injections**  
In virtue of all the advantages listed here "MN\*" needles revolutionize our aesthetic approach to the face and transform the minds of our patients who in the future will consult us with greater serenity.

# Cannula

Main advantages & skin tissue



## Anteis cannula range

| Ref. Cannulas    | Anteis Gel  |  Manual Injection |  AIS Injection + RTU** |
|------------------|---|---|---|
| PNE 010030 (MN*) | Esthélis (Basic, Soft),<br>Jolidermis 18<br>Jolidermis 24 | 27G * 37 mm   | 27G * 37 mm or<br>30G * 27 mm   |
| PNE 010010 (MN*) | Fortélis, Jolidermis 24+,<br>Jolidermis 25                | 25G * 40 mm   | 25G * 40 mm or<br>27G * 37 mm   |
| PNE 020010       | Modélis, Jolidermis Volume                                | 22 G * 50 mm  | 25G * 40 mm   |

Thinner needles thanks to AIS

\* MN: Magic Needle™ · \*\* RTU: Ready-to-use syringes (adapted at origin for AIS)

## Ideal injection setup

Combine Anteis Injection System with cannula and increase the comfort of your patients